

Komet® ZR Flash Polishers™

For Dental Labs



The Best Solution For High-Strength Ceramics

Komet ZR Flash Polishers have been developed to efficiently impart a high shine when used for the final polishing of restorations fabricated of zirconium oxide and other hard, high-strength ceramic materials, including lithium disilicate, aluminum oxide, and pressed ceramics. ZR Flash Polishers incorporate interspersed diamonds for consistently reliable high performance.

The two-step polishers are identified by color: blue for smoothing and pre-polishing and light gray for high-shine final polishing, each available in four sizes and configurations: a knife-edged 26-mm-diameter version and flatedge 17-mm and 11-mm-diameter versions.

The laboratory line of polishers for high-strength ceramics now includes a 16.5-mm-long, tapered flame-shape ZR Flash Polisher, which also is offered in blue (for smoothing and pre-polishing) and light-gray (for final polishing) versions.

Pre-polishing

94011C.HP.260

94012C.HP.110

94013C.HP.170

94018C.HP.055

94019C.HP.110

High-shine polishing


94011F.HP.260

94012F.HP.110

94013F.HP.170

94018F.HP.055

94019F.HP.110



ZR Flash Polishers For Chairside

To meet the unique preliminary and final polishing requirements of high-strength ceramic materials at chairside, Komet ZR Flash Polishers have been developed to smooth, polish, and impart a high shine when applied to restorations fabricated of zirconium oxide, aluminum oxide, lithium disilicate, and pressed ceramics.

Incorporating interspersed diamonds for consistently reliable performance, the polishers produce outstanding results in two steps, each of which can be easily identified by the color of the Komet ZR Flash Polisher: blue for smoothing and pre-polishing and light gray for high-shine final polishing. Available in styles, sizes, and shapes designed specifically for intraoral applications, Komet ZR Flash Polishers for chairside use include two flame-shaped versions (each available in blue and light gray) as well as a cup configuration (also offered in blue and light-gray versions).

A kit (LD0707A) containing two flame-shaped and two cup-configured ZR Flash Polishers as well as three ZR-Diamonds™ – all the elements for effective chairside polishing and adjustment of high-strength ceramic restorations – also is available from Komet USA.



94020C.RA.040




94020F.RA.040



94021C.RA.050



94021F.RA.050



94022C.RA.060



94022F.RA.060



Kit Order LD0707A

For questions or
to place an order



888-KOMETUSA



info@kometusa.com



kometusa.com

